

FKBP12 Ligand-Linker Conjugate 2

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	FKBP12 Ligand-Linker Conjugate 2 is a conjugate composed of an FKBP12 target protein ligand and a linker. It is utilized in the synthesis of FKBP12PROTACFM4.
Targets(IC50)	Target Protein Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481