

## GLP-1R agonist 29

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	GLP-1R agonist 29 (Compound 20) is a GLP-1R agonist that induces hGLP-1R-mediated cAMP stimulation with an EC50 of 0.018 nM. It exhibits favorable pharmacokinetic properties and shows good in vivo exposure, with an AUC <sub>0-∞,sc</sub> of 77688 ng·h/mL.
Targets(IC50)	Glucagon Receptor

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481