

EGFR-IN-153

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	EGFR-IN-153 (Compound 5l) is an EGFR inhibitor that suppresses the proliferation of MDA-MB-231 cells with an IC50 of 0.73 $\mu$ M and induces apoptosis. It is applicable for research in breast cancer.
Targets(IC50)	Apoptosis,EGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481