

pGk13a

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	pGk13a is an amphiphilic membrane labeling probe containing azido groups (which can bind with fluorophores). It enables high-resolution imaging of cell membranes in ultra-micro expansion microscopy (umExM), aiding in the observation of membrane-associated structures and proteins. pGk13a is applicable in studies of neuronal structures.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481