

XW-032

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	XW-032 is an apo-IDO1 inhibitor with an IC50 of 21 nM. In the CT26 syngeneic mouse model, XW-032 demonstrates significant in vivo antitumor efficacy with a TGI of 63%, showing potential for cancer research applications.
Targets(IC50)	IDO,Indoleamine 2,3-Dioxygenase (IDO)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481