

I-67

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	I-67, a structural isomer of I-66, is a potent β -catenin inhibitor with an IC50 of less than 50 nM. It disrupts the interaction between β -catenin and T cell factor (TCF) transcription factors and exhibits anticancer properties.
Targets(IC50)	Wnt/beta-catenin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481