

DL-01 TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DL-01 TFA is the trifluoroacetate salt form of DL-01. It is a drug-linker conjugate designed for use in the synthesis of antibody-drug conjugates (ADC).
Targets(IC50)	Drug-Linker Conjugates for ADC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481