

## Carbonic anhydrase inhibitor 29

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Carbonic anhydrase inhibitor29 (Compound 5d) is an inhibitor targeting carbonic anhydrase IX and XII, with an inhibition constant (Ki) of 26.6 nM for carbonic anhydrase IX and a Ki of 10.9 nM for carbonic anhydrase XII. Carbonic anhydrase inhibitor29 can be used in cancer research.
Targets(IC50)	Carbonic Anhydrase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481