

Top1-IN-2

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Top1-IN-2 (Compound 1a) is an inhibitor of topoisomerase 1 (Top1). It suppresses the growth of P-gp drug-resistant tumor cells and induces apoptosis.
Targets(IC50)	Apoptosis,Topoisomerase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481