

## VH 101-amide-piperidine-Pip-alkyne

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	VH 101-amide-piperidine-Pip-alkyne is an E3LigaseLigand-Linker Conjugate used in the synthesis of PROTAC degraders. The linker component is tert-Butyl 4-ethynyl-[1,4'-bipiperidine]-1'-carboxylate, and the VHL-based E3 ubiquitin ligase ligand is VH 101, acid.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481