

FITC-PEG-MAL (MW 1000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | FITC-PEG-MAL (MW 1000) is a fluorescent dye composed of polyethylene glycol (PEG), FITC, and Maleimide (MAL). It is applicable in fluorescence labeling and imaging. |
| Targets(IC50) | Others   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481