

Bemiltenase alfa

Chemical Properties

CAS No. : 2766038-73-5

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Bemiltenase alfa is a hemostatic agent that targets coagulation factor X (FX). It activates FX to facilitate the conversion of prothrombin into thrombin, which subsequently transforms fibrinogen into fibrin, thereby forming a stable blood clot and exerting its hemostatic effect. Bemiltenase alfa shows potential for research in bleeding symptoms associated with hemophilia.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481