

## (R)-PIN1 ligand-2

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Powder: -20°C for 3 years | In solvent: -80°C for 1 year**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	(R)-PIN1 ligand-2 is a PIN1 ligand that can be utilized in the synthesis of PROTACs, such as P1D-34.
Targets(IC50)	Ligands for Target Protein for PROTAC

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481