

AS-0017445

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | AS-0017445 is an inhibitor targeting the main protease of the coronavirus and the virus causing Middle East Respiratory Syndrome (MERS). It functions by inhibiting viral protein processing within host cells, thereby preventing viral replication. |
| Targets(IC50) | SARS-CoV, Virus Protease |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481