

PROTAC AR Degradar-10

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	PROTACAR Degradar-10 is a protein degrader targeting androgen receptors (AR) with a DC50 value of ≤ 100 nM, and can be used in prostate cancer research.
Targets(IC50)	Androgen Receptor,PROTACs

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481