

QI-18

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	QI-18 is a fibroblast activation protein (FAP) inhibitor with an IC50 value of 0.50 nM, offering a 6.5-fold increase in activity compared to UAMC-1110 (IC50 3.25 nM). It can be utilized in the synthesis of highly tumor-selective and high-dose radiotracers for the diagnosis and treatment of tumors.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481