

DAO-dBET1

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DAO-dBET1, a dual-action PROTAC incorporating the PROTAC degrader dBET1, effectively serves as a potent BRD4 degrader and exhibits a DC50 value of 277 nM in the presence of CoCl <sub>2</sub> . Additionally, DAO-dBET1 inhibits hypoxia and Cath-L activation with an IC50 value of 281 nM.
Targets(IC50)	Epigenetic Reader Domain,PROTACs

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481