

DSPE-PEG2000-LTLRWVGLMS

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | DSPE-PEG2000-LTLRWVGLMS is a PEG compound composed of DSPE and the decapeptide (LTLRWVGLMS). Chondroitin sulfate proteoglycan NG2 is the receptor for LLRWVGLMS. LLRWVGLMS demonstrates homing to pericytes associated with tumor vasculature. DSPE-PEG2000-LTLRWVGLMS can be utilized for drug delivery. |
| Targets(IC50) | Liposome |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481