

FITC-PEG-FA (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | FITC-PEG-FA (MW 1000) is a fluorescent dye composed of folic acid (FA), polyethylene glycol (PEG), and FITC. Through specific chemical reactions, fluorescein is linked with PEG and folic acid to create a molecule with unique properties. FITC-PEG-FA (MW 1000) is suitable for fluorescent labeling and imaging. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481