

## Armenian hamster IgG1 kappa, Isotype Control

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Armenian hamster IgG1 kappa, Isotype Control, is a monoclonal antibody derived from the Armenian hamster, serving as the isotype control for Armenian hamster IgG1κ antibodies.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481