

(R,S,S)-VH032-Me-N-Boc-7-aminoheptanoic acid

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	(R,S,S)-VH032-Me-N-Boc-7-aminoheptanoic acid is a conjugate of an E3 ligase ligand and a linker (E3LigaseLigand-Linker Conjugates), composed of N-Boc-7-aminoheptanoic acid and the corresponding linker. It is utilized in the synthesis of dTAGV-1-NEG.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481