

FITC-PEG-COOH (MW 400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	FITC-PEG-COOH (MW 400) is a fluorescent dye comprised of polyethylene glycol (PEG), FITC, and a carboxyl group (COOH). FITC-PEG-COOH (MW400) is suitable for fluorescent labeling and imaging.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481