

## Kv7.2/Kv7.3 modulator-2

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Kv7.2/Kv7.3 modulator-2 (compound 7d) is a Kv7.2/Kv7.3 channel modulator with a pEC50 value of 7.42 for activating the Kv7.2/3 channels. This compound is applicable in the research of epilepsy, depression, brain injury, and pain.
Targets(IC50)	Potassium Channel

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481