

## Heparin-Cy3

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Powder: -20°C for 3 years | In solvent: -80°C for 1 year**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Heparin-Cy3 is a fluorescent dye composed of Cy3 and Heparin. It is used for cell tracking and drug delivery.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481