

## DOPE-DOTA sodium

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

DOPE-DOTA sodium is a chelating lipid that serves as a crucial contrast agent in magnetic resonance imaging (MRI). It is applicable in studies related to cancer, blood-brain barrier disruption, aneurysms or plaque buildup, arthritis inflammation, as well as liver function and lesions.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481