

DSPE-PEG1000-CTT2

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG1000-CTT2 is a PEG compound comprised of DSPE and a gelatinase inhibitor [CTT2 (CTTHWGFTLC)]. The CTT2 (CTTHWGFTLC) peptide is capable of specifically targeting tumors. DSPE-PEG1000-CTT2 can be utilized for drug delivery.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481