

FITC-PEG-FA (MW 10000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	FITC-PEG-FA (MW 10000) is a fluorescent dye composed of Folic acid (FA), polyethylene glycol (PEG), and FITC. Through specific chemical reactions, fluorescein is conjugated with PEG and Folic acid, resulting in a molecule with unique properties. FITC-PEG-FA (MW 10000) is suitable for fluorescence labeling and imaging.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481