

DSPE-PEG5000-TH

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG5000-TH is a PEG compound formed from DSPE and the pH-responsive cell-penetrating peptide (TH). In acidic environments like the tumor microenvironment, TH is activated, allowing it to selectively deliver small molecules, oligonucleotides, proteins, and similar substances into tumor cells.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481