

HRG038

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HRG038 is a selective covalent BRD4 degrader. It binds covalently to the E3 ubiquitin ligase CUL4 DCAF16, recruiting it to BRD4, which results in BRD4 ubiquitination and subsequent proteasomal degradation. HRG038 holds potential for cancer research.
Targets(IC50)	Epigenetic Reader Domain

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481