

BRD4 ligand 7

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	BRD4 ligand 6 is a PROTAC target protein ligand (Ligand for Target Protein for PROTAC). It serves as a key component in synthesizing the PROTAC CBP/BRD4 degrader-1.
Targets(IC50)	Ligands for Target Protein for PROTAC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481