

DSPE-PEG5000-A7R

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG5000-A7R is a PEG compound composed of DSPE and the tumor vascular targeting peptide (A7R). The peptide A7R exhibits high affinity and specificity for VEGFR-2, which is overexpressed in various tumors.
Targets(IC50)	Liposome, VEGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481