

C8 DG PEG 2000

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

C8 DG PEG 2000 (1,2-dioctanoyl-rac-glycero-3-methoxypolyethylene glycol-2000) is a PEG lipid consisting of a glycerol backbone with two C8 acyl chains and a PEG 2000 unit. It is applicable in the formation of lipid nanoparticles and liposomes.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481