

m-PEG1000-PEI-Cy5

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	mPEG-PEI-Cy5 (MW 1000) is a fluorescent dye composed of methoxy polyethylene glycol (mPEG), Cy5, and polyethyleneimine (PEI). This compound is suitable for applications in fluorescence labeling and imaging.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481