

DSPE-PEG1000-NGR

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG1000-NGR is a PEG compound composed of DSPE and the Asn-Gly-Arg (NGR) peptide. The NGR peptide has the capability to target the CD13 antigen present in tumor vasculature. DSPE-PEG1000-NGR can be utilized for drug delivery.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481