

## PROTAC STING degrader-3

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PROTACSTING degrader-3 (Compound ST9) is a STINGPROTAC-type degrader with a DC50 of 0.62 $\mu$ M. It induces STING degradation via the ubiquitin-proteasome pathway. This compound exerts anti-inflammatory effects by inhibiting the STING/TBK1/NF- $\kappa$ B signaling. Additionally, it offers renal protection and can be used in research on acute kidney injury.
Targets(IC50)	NF- $\kappa$ B,PROTACs,STING

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481