

Phytolacca americana Lectin

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Phytolacca americana Lectin (PWM) is a lectin specific for carbohydrates containing N-acetylglucosamine. It stimulates peripheral lymphocytes to undergo mitosis by binding to the cell surface. Phytolacca americana Lectin can serve as a probe for specific binding to biomolecules and is useful as a biomaterial or organic compound in life science research.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481