

Phaseolus vulgaris Lectin-P

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Phaseolus vulgaris Lectin-P is a protein form of Phytohemagglutinin obtained prior to isolating and purifying Phytohemagglutinin-L (PHA-L) and PHA-E. It can be used to stimulate the proliferation of peripheral mononuclear cells, enhance the production of certain cytokines, and promote the expression of membrane surface proteins.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481