

dGEM3

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	dGEM3 is a molecular glue degrader targeting GEMIN3 (a component of the survival of motor neuron complex). It is applicable in research on diseases related to abnormal protein regulation, such as cancer and neurodegenerative disorders.
Targets(IC50)	Molecular Glues,PROTACs

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481