

Mc-PEG4-Val-Ala-PAB-Exatecan

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Mc-PEG4-Val-Ala-PAB-Exatecan is a drug-linker conjugate used in ADC (antibody-drug conjugate) formulations. This compound comprises an ADC linker component (Mc-PEG4-Val-Ala-PAB) and the DNA topoisomerase I inhibitor, Exatecan. |
| Targets(IC50) | Drug-Linker Conjugates for ADC,Topoisomerase |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481