

DSPE-PEG2000-pPB

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG2000-pPB is a PEG compound comprised of DSPE and the cyclic oligopeptide (pPB). The pPB component exhibits a high binding affinity for PDGFR $\beta$ , which is overexpressed on activated hepatic stellate cells (HSC). DSPE-PEG2000-pPB is utilizable for drug delivery applications.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481