

DSPE-PEG2000-K237

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG2000-K237 is a PEG compound consisting of DSPE and the KDR-targeting peptide (K237). K237 functionally disrupts the interaction between VEGF and KDR receptors, resulting in significant biological effects such as inhibition of angiogenesis and tumor growth. DSPE-PEG2000-K237 is applicable in drug delivery.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481