

DSPE-PEG5000-ANG

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG5000-ANG is a PEGylated compound formed by DSPE and the dual-targeting ligand (Angiopep-2, ANG). ANG exhibits high binding efficacy to LRP1 and has been utilized for targeting glioma delivery. DSPE-PEG5000-ANG is applicable in drug delivery.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481