

## Cholesteryl-PEG1000-Biotin

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Cholesteryl-PEG1000-Biotin is a lipid product with a cholesterol backbone linked to a PEG unit on one end and a biotin unit on the other. Biotin (Biotin) is an enzyme cofactor useful for protein labeling, while PEG is a low-toxicity hydrophilic polymer. Cholesteryl-PEG1000-Biotin is applicable in drug circulation time studies.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481