

DSPE-PEG5000-MPG

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG5000-MPG is a PEG compound composed of DSPE and a peptide carrier (MPG). MPG is derived from the nuclear localization sequence (NLS) of the SV40 large T antigen and the fusion peptide domain of HIV-1 gp41. This compound effectively delivers short oligonucleotides into cells independently of the endosomal pathway.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481