

XZ338

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	XZ338 is a highly selective degrader targeting BCL-XL without degrading BCL-2. It exhibits an IC50 value of 3.7 nM against MOLT-4 cells and possesses antiproliferative properties, making it useful for cancer research.
Targets(IC50)	Bcl-2 Family

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481