

DSPE-PEG1000-M2pep

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG1000-M2pep is a PEG compound comprised of DSPE and the M2 macrophage-targeting peptide (M2pep). M2pep acts as a mono- or multivalent ligand conjugated with pro-apoptotic peptides to focus the peptides' toxic effects on M2 macrophages. DSPE-PEG1000-M2pep can be utilized for drug delivery.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481