

DSPE-PEG2000-R6H4

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG2000-R6H4 is a PEG compound composed of DSPE and a pH-responsive cell-penetrating peptide (R6H4). The R6H4 component is utilized for pH-responsive anticancer drug delivery. Additionally, DSPE-PEG2000-R6H4 is applicable in drug delivery systems.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481