

Rapavir

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Rapavir (JH-B10) is a selective and potent inhibitor of sodium taurocholate co-transporting polypeptide (NTCP), with an IC50 value of 1.8 nM for taurocholic acid-d4 (TCA-d4) uptake. It exhibits antiviral activity against HBV infection by directly binding to NTCP and preventing the virus from entering cells. Rapavir is promising for HBV infection research.
Targets(IC50)	HBV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481