

BDM_92499

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BDM_92499 is a nanomolar, selective IRAP inhibitor with an IC50 of 3.4 nM. It also inhibits ERAP1 and ERAP2, with IC50 values of 0.46 μM and 4.2 μM, respectively.
Targets(IC50)	Aminopeptidase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481