

## PROTAC KDM4 degrader-1

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PROTACKDM4 degrader-1 (Compound 11) is a proteolysis-targeting chimera (PROTAC) degrader specifically designed for KDM4. It effectively degrades KDM4A-C while sparing KDM4D. In esophageal cancer cells, PROTACKDM4 degrader-1 demonstrates strong antiproliferative activity, inducing apoptosis and cell cycle arrest, and inhibits lysine demethylation of histone H3.
Targets(IC50)	Apoptosis,Histone Demethylase,PROTACs

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481